

G1061 pFBAAVTREmiSafeBgHpA CMVrtTAwtIRESegFP SV40pA Inducible miRNA Plasmid Control



Plasmid Features:

<u>Coordinates</u>	<u>Feature</u>
194-348	Tn7L
376-511	SV40pA (complementary -outside ITRs)
678-818	AAV2 ITR (141bp)
878-1193	TRE Tight promoter Clonotech
1211-1296	miSafe
1325-1556	BgHpA
1632-2211	CMV Promoter
2240-2992	rtTA Advanced Tet On Transactivator Clonotech
3007-3593	wtIRES
3594-4313	EGFP
4386-4522	SV40pA
4594-4723	AAV2 ITR (130bp)
5280-5813	Gentamicin
5880-6104	Tn7R
7168-8028	B-lactamase (Ampicillin) Complementary


Antibiotic Resistance: Ampicillin and Gentamicin
 Bacterial Backbone: pFastBac™ Invitrogen
 TREtight: Clonotech
 rtTA Advanced: Clonotech

Note: To check the integrity of the AAV ITR's perform single restriction enzyme digestions with XmaI, BssHII, and MscI.

This vector was designed with EGFP under the constitutive expression from the CMV promoter for tracking while the expression of the miRNA cassette can be inducibly driven by the tetracycline-regulated TRE promoter for RNA interference (RNAi) analysis in mammalian cells.

miSafe is a scrambled control developed by Ryan Boudreau in the lab of Beverly Davidson.

Please reference: Boudreau, R. L., et al. (2011). "Rational Design of Therapeutic siRNAs: Minimizing Off-targeting Potential to Improve the Safety of RNAi Therapy for Huntington's Disease." Mol Ther 19(12): 2169-2177.




>G1061 pFBAAVTREmiSafeBgHpACMvrtAwtIRESegFPSV40pA

CCATTCGCCATTCAGGCTGCAAAATAAGCGTTGATATTCAGTCAATTACAAACATTAATAACGAAGAGATGACAGAAAA
ATTTTCATTCTGTGACAGAGAAAAAGTAGCCGAAGATGACGGTTTGTACATGGAGTTGGCAGGATGTTTGATTA
AACATAACAGGAAGAAAAATGCCCGCTGTGGGCGGACAAAATAGTTGGGAAGTGGGAGGGGTGGAAATGGAGTTTT
TAAGGATTATTTAGGGAAGAGTGACAAAATAGATGGGAAGTGGGTGTAGCGTCGTAAGCTAATACGAAAATTA
TGACAAAATAGTTTGGAACTAGATTTCACTTATCTGGTTCCGATCTCCTAGGCGATATCAGTGATCAGATCCAGACAT
GATAAGATACATTGATGAGTTTGGACAAACCACAAGTGAATGCAGTGAAAAAATGCTTTATTTGTGAAATTTGTG
ATGCTATTGCTTTATTTGTAACCATTATAAGCTGCAATAAACAAGTTAACACAACAATTGCATTCATTTTATGTTTC
AGGTTACAGGGGGAGGTGTGGGAGGTTTTTTAAAGCAAGTAAAACCTCTACAAATGTGGTATGGCTGATTATGATCCTC
TAGTACTTCTCGACAAGCTTACATTATTGAAGCATTTATCAGGGTTATTGTCTCAGACTGCGAGGCTGCGCGCTCG
CTCGTCACTGAGGCCGCCGGGCAAAGCCCGGGCGTCGGGCGACTTTGGTGGCCGGCTCAGTGAGCGAGCGAGCGC
GCAGAGAGGGAGTGGCCAATCCATCACTAGGGTTCTTGTAGTTAATGATTAACCCGCCATGCTACTTATCTACGTA
GCCATGCTCTAGTGAATTCGAGTTTACTCCCTATCAGTGATAGAGAACGTATGTGAGTTTACTCCCTATCAGTGATA
GAGAACGATGTGAGTTTACTCCCTATCAGTGATAGAGAACGTATGTGAGTTTACTCCCTATCAGTGATAGAGAACG
TATGTGAGTTTACTCCCTATCAGTGATAGAGAACGTATGTGAGTTTACTCCCTATCAGTGATAGAGAACGTATGTG
AGTTTACTCCCTATCAGTGATAGAGAACGTATGTGAGTTTACTCCCTATCAGTGATAGAGAACGTATGTGAGTTT
TTTACTCCCTATCAGTGATAGAGAACGTATGTGAGTTTACTCCCTATCAGTGATAGAGAACGTATGTGAGTTTACT
GTAAGCCACAGATGGGTAACCGCTAAGTCGTTTCGCTACGCCCTACTAGAAGCGAGCGCCACCGCGGACAGCGGGGAGATCC
AGACATGATAAGATCCTCGACTGTGCCTTCTAGTTGCCAGCCATCTGTTGTTTGGCCCTCCCCGTCCTTCTTGACCC
TGGAAGGTGCCACTCCCCTGTCCTTCCCTAATAAAATGAGGAAATTGCATCGCATTGCTGAGTAGGTGTCATTCTA
TTCTGGGGGGTGGGGTGGGGCAGGACAGCAAGGGGGAGGATTGGGAAGACAATAGCAGGCATGCTGGGGAAGCCGGTC
GCAGGCGTAGGGTAGGTTTGGAGCAGATAAACCATTAAGTGATTAACAACGACGGCAGGACGCGGCCGTTATTAATAGT
AATCAATTACGGGGTCATTAGTTCATAGCCCATATATGGAGTTCGCGTTACATAACTTACGGTAAATGGCCCGCCTG
GCTGACCGCCCAACGACCCCCGCCATTGACGTCAATAATGACGTATGTTCCCATAGTAACGCCAATAGGGACTTTCCA
TTGACGTCAATGGGTGGAGTATTTACGGTAAACTGCCACTTGGCAGTACATCAAGTGTATCATATGCCAAGTACGCC
CCCTATTGACGTCAATGACGGTAAATGGCCCGCCTGGCATTATGCCAGTACATGACCTTATGGGACTTTCCCTACTTGG
CAGTACATCTACGTATTAGTCATCGCTATTACCATGGTGATGCGGTTTTGGCAGTACATCAATGGGCGTGGATAGCGG
TTTGACTCACGGGGATTTCGAAGTCTCCACCCCATGACGTCAATGGGAGTTTGTGTTTGGCACCAAAATCAACGGGACT
TTCCAAAATGTGTAACAACCCGCCCATGACGCAAAATGGGCGGTAGGCGTGTACGGTGGGAGGTCTATATAAGCA
GAGCTGGTTTGTGAAACCGTCAGATCCGCTAGCTACCGGTTAATTAAGCCACCATGTCCAGACTGGACAAGAGCAAAG
TCATAAACGGCGCTCTGGAATTACTCAATGGAGTCGGTATCGAAGGCCTGACGACAAGGAAACTCGCTCAAAGCTGG
GAGTTGAGCAGCCTACCCTGTACTGGCACGTGAAGAACAAGCGGGCCCTGCTCGATGCCCTGCCAATCGAGATGCTGGA
CAGGCATCATACCCACTTCTGCCCCCTGGAAGGCGAGTCATGGCAAGACTTCTGCGGAACAACGCCAAGTCATTCCGC
TGTGCTCTCCTCTCACATCGCGACGGGGCTAAAGTGCATCTCGGCACCCGCCAACAGAGAAACAGTACGAAACCTGG
AAAATCAGCTCGCGTTCTGTGTCAGCAAGGCTTCTCCCTGGAGAACGCACTGTACGCTCTGTCCGCCGTGGGCCACTT
TACACTGGGCTGCGTATTGGAGGAACAGGAGCATCAAGTAGCAAAAGAGGAAAGAGAGACACCTACCACCGATTCTAT
GCCCCACTTCTGAGACAAGCAATTGAGCTGTTGACCGGCAGGGAGCCGAACCTGCCTTCTTTTGGCCTGGAACATA
ATCATATGTGGCCTGGAGAAACAGCTAAAGTGCAGAAAGCGGCGGGCCGGCGACGCCCTTGACGATTTTGACTTAGAC
ATGCTCCCAGCCGATGCCCTTGACGACTTTGACCTTGATATGCTGCTGCTGACGCTTTGACGATTTTGACCTTGACA
TGCTCCCCGGGTAAGTCGACTTGAATTCGCCCTCTCCCTCCCCCCCCCTAACGTTACTGGCCGAAGCCGCTTGAATA
AGGCCGGTGTGCGTTTGTCTATATGTTATTTCCACCATATTGCCGTCTTTTGGCAATGTGAGGGCCCGAAACCTGGC
CCTGTCTTCTTGACGAGCATTCTAGGGGTCTTTCCCTCTCGCAAAGGAATGCAAGGTCTGTTGAATGTCGTGAAGG
AAGCAGTTCCTCTGGAAGCTTCTGAAAGACAAAACAGTCTGTAGCGACCCTTTGCAGGCAGCGGAACCCCCACCTGG
CGACAGGTGCCTCTGCGGCCAAAAGCCACGTGTATAAGATACACCTGCAAAGGCGGCACAACCCAGTGCCACGTTGTG
AGTTGGATAGTTGTGAAAGAGTCAAAATGGCTCTCCTCAAGCGTATTCAACAAGGGGCTGAAGGATGCCCAGAAGGTA
CCCCATTGTATGGGATCTGATCTGGGGCTCGGTACACATGCTTTACATGTGTTTGTGAGGTTAAAAACGTCTAG
GCCCCCGAACCACGGGGACGTGGTTTTCTTTGAAAAACACGATGATAATATGGCCACAACCATGGTGAGCAAGGGCG
AGGAGCTGTTACCGGGGTGGTGGCCATCCTGGTGCAGCTGGACGGCGACGTAACGGCCACAAGTTCAGCGTGTCCGG
CGAGGGCGAGGGCGATGCCACCTACGGCAAGCTGACCTGAAGTTCATCTGCACCACCGGCAAGCTGCCCGTGGCCCTGGC
CCACCTCGTGACCACCTGACCTACGGCGTGCAGTGCTTCAGCCGCTACCCCGACCACATGAAGCAGCAGGACTTCTTC
AAGTCCGCCATGCCGAAGGCTACGTCCAGGAGCGCACCATCTTCTTCAAGGACGACGGCAACTACAAGACCCCGCGCGA
GGTGAAGTTGAGGGCGACACCTGGTGAACCGCATCGAGCTGAAGGGCATCGACTTCAAGGAGGACGGCAACATCCTG
GGGCACAAGCTGGAGTACAACATAACAGCCACAACGTCTATATCATGGCCGACAAGCAGAAGAACGGCATCAAGGTG
AACTTCAAGATCCGCCACAACATCGAGGACGGCAGCGTGCAGCTCGCCGACCCTACCAGCAGAACACCCCATCGGCGA
CGCCCCGTGCTGCTGCCGACAACCACTACCTGAGCACCCAGTCCGCCCTGAGCAAAGACCCCAACGAGAAGCGCGATC
ACATGGTCTGCTGGAGTTCGTGACCGCCCGGGATCACTCTCGGCATGGACGAGCTGTACAAGTAAATCGATTTACG
CGTGAGCTCCAATTCGCCCTATAGTGAGTCGATTACGCGCGCAGCGGCCGACCATGGCCCAACTTGTATTGACAGT

TATAATGGTTACAAATAAAGCAATAGCATCACAAATTTACAAATAAAGCATTTTTTTTCACTGCATTCTAGTTGTGGT
TTGTCCAAACTCATCAATGTATCTTATCATGTCTGGATCTCCGGACACGTGCGGACCGAACTAGAGCATGGCTACGTAG
ATAAGTAGCATGGCGGGTTAATCATTAACACAAAGGAACCCCTAGTGATGGAGTTGGCCACTCCCTCTCTGCGCGCTCG
CTCGCTCACTGAGGCCGGGGACCAAAGGTCGCCGACGCCGGGGCGCCTCAGTGAGCGAGCGAGCGCGCAGCTGCCTG
CAGGCATGCAAGCTGTAGCCAACCACTAGAACTATAGCTAGAGTCTGGGCGAACAAACGATGCTCGCCTTCCAGAAAA
CCGAGGATGCGAACCCTTCCATCCGGGGTCAGCACCACCGGCAAGCGCCGACGCGCCGAGGTCTTCCGATCTCCTGAAG
CCAGGGCAGATCCGTGCACAGCACCTTGCCGTAGAAGAACAGCAAGGCCGCAATGCCTGACGATGCGTGGAGACCGAA
ACCTTGCCTCGTTCGCCAGCCAGGACAGAAATGCCTCGACTTCGCTGCTGCCAAGGTTGCCGGGTGACGCACACCGTG
GAAACGGATGAAGGCACGAACCCAGTTGACATAAGCCTGTTCCGGTTCGTAAACTGTAATGCAAGTAGCGTATGCGCTC
ACGCAACTGGTCCAGAACCTTGACCGAACGCAGCGGTGGTAACGGCGCAGTGCGGTTTTTCATGGCTTGTATGACTGT
TTTTTTGTACAGTCTATGCCTCGGGCATCCAAGCAGCAAGCGGTTACGCCGTGGGTGATGTTTGTATGTTATGGAGC
AGCAACGATGTTACGCAGCAGCAACGATGTTACGCAGCAGGGCAGTCGCCCTAAAACAAAGTTAGGTGGCTCAAGTAT
GGGCATCATTGACATGTAGGCTCGGCCCTGACCAAGTCAAATCCATGCGGGCTGCTCTTGATCTTTTCGGTCTGAG
TTCGGAGACGTAGCCACCTACTCCAACATCAGCCGACTCCGATTACCTCGGGAACCTGCTCCGTAGTAAGACATTCA
TCGCGCTTGCTGCCTTCGACCAAGAAGCGGTTGTTGGCGCTCTCGCGCTTACGTTCTGCCCAAGTTTGAGCAGCCGG
TAGTGAGATCTATATCTATGATCTCGCAGTCTCCGGCGAGCACCGGAGGCAGGGCATTGCCACCGCGCTCATCAATCTC
CTCAAGCATGAGGCCAACCGCTTGGTGTATGTGATCTACGTGCAAGCAGATTACGGTGACGATCCCGCAGTGGCTC
TCTATACAAAGTTGGGCATACGGGAAGAAGTGATGCACTTTGATATCGACCAAGTACCGCCACCTAACAAATTCGTTCA
AGCCGAGATCGGCTTCCCGCCCGCGGAGTTGTTCCGTAATTTGTCACAACGCCGGAATATAGTCTTTACCATGCCCTT
GGCCACGCCCTCTTAATACGACGGCAATTTGCACTTCAGAAAATGAAGAGTTTGTCTTAGCCATAACAAAAGTCCA
GTATGCTTTTTACAGCATAACTGGACTGATTTCAAGTTACAACCTATTCTGTCTAGTTTAAAGACTTTATTGTCAAGT
TTAGATCTATTTTGTTCAGTTTAAAGACTTTATTGTCCGCCACACCCGCTTACGCAGGGCATCCATTTATTACTCAACC
GTAACCGATTTTGCAGGTTACGCGGCTGGTCTATGCGGTGTGAAATACCGCACAGATGCGTAAGGAGAAAAATACCGC
ATCAGGCGCTCTTCCGCTTCCCTCGCTCACTGACTCGCTGCGCTCGGTCTCGGCTGCGGCGAGCGGTATCAGCTCACTC
AAAGGCGGTAATACGGTTATCCACAGAATCAGGGGATAACGCAGGAAAAGAAATGTGAGCAAAAAGGCCAGCAAAAAGGC
CAGGAACCGTAAAAAGGCCGCGTTGCTGGCGTTTTTCCATAGGCTCCGCCCCCTGACGAGCATCACAAAAATCGACGC
TCAAGTCAGAGGTGGCGAAACCCGACAGGACTATAAAGATACCAGGCGTTTTCCCTGGAAGCTCCCTCGTGCGCTCTC
CTGTTCCGACCCTGCCGCTTACCGGATACCTGTCCGCTTTCTCCCTTCGGGAAGCGTGGCGCTTTCTCATAGCTCACGC
TGTAGGTATCTCAGTTCGGTGTAGGTGCTTCCGCTCCAAGCTGGGCTGTGTGCACGAACCCCCGTTACGCCCCACCGCT
GCGCCTTATCCGGTAACTATCGTCTTGAGTCCAACCCGTAAGACACGACTTATCGCCACTGGCAGCAGCCACTGGTAA
CAGGATTAGCAGAGCGAGGTATGTAGGCGGTGCTACAGAGTTCCTGAAAGTGGTGGCCTAACTACGGCTACACTAGAAG
AACAGTATTTGGTATCTGCGCTCTGCTGAAGCCAGTTACCTTCGGAAAAAGAGTTGGTAGCTCTTGATCCGGCAAAACA
AACCACCGCTGGTAGCGGTGGTTTTTTTTGTTTGAAGCAGCAGATTACGCGCAGAAAAAAGGATCTCAAGAAGATCC
TTTGATCTTTTCTACGGGGTCTGACGCTCAGTGAACGAAAACCTCACGTTAAGGGATTTTGGTCATGAGATTATCAA
AAGGATCTTACCTAGATCTTTTAAATTAATAAATGAAGTTTTAAATCAATCTAAAGTATATATGAGTAAACTTGGTC
TGACAGTTACCAATGCTTAATCAGTGAGGCACCTATCTCAGCGATCTGTCTATTTTCGTTTATCCATAGTTGCCTGACTC
CCCCTCGTGTAGATAACTACGATACGGGAGGGCTTACCATCTGGCCCCAGTGTGCAATGATACCGCGAGACCCACGCT
CACCGCTCCAGATTTATCAGCAATAAACCCAGCCAGCCGGAAGGGCCGAGCGCAGAAGTGGTCTGCAACTTTATCCGC
CTCCATCCAGTCTATTAATTTGTTGCCGGGAAGCTAGAGTAAGTAGTTCCGCAAGTTAATAGTTTGCGAACGTTGTTGC
CATTGCTACAGGCATCGTGGTGTACGCTCGTCTTGGTATGGCTTATTACGCTCCGTTCCCAACGATCAAGGCGA
GTTACATGATCCCCATGTTGTGCAAAAAGCGGTTAGCTCCTTCGGTCTCCGATCGTTGTGAGAAGTAAGTTGGCCG
CAGTGTATCACTCATGGTTATGGCAGCACTGCATAATCTCTTACTGTGATGCCATCCGTAAGATGCTTTTCTGTGAC
TGGTGAGTACTCAACCAAGTCAATCTGAGAATAGTGTATGCGGCGACCGAGTTGCTCTTGGCCGGCGTCAATACGGGA
TAATACCGCGCCACATAGCAGAACTTTAAAAGTGCTCATCATTGGAAAACGTTCTTCCGGGGCGAAAACCTCTCAAGGAT
CTTACCGCTGTTGAGATCCAGTTCGATGTAACCCACTCGTGACCCAACTGATCTTACGATCTTTTACTTTTACCAGC
GTTTCTGGGTGAGCAAAAACAGGAAGGCAAAATGCCGCAAAAAGGGAATAAGGGCGACACGGAAATGTTGAATACTC
ATACTCTTCTTTTTCAATATTTATGAAGCATTATCAGGGTTATTGTCTCATGAGCGGATACATATTTGAATGTATT
TAGAAAAATAACAAATAGGGGTTCCGCGCACATTTCCCGGAAAAGTGCCACCTAAATGTAAGCGTTAATATTTTGT
TAAAATTCGCTTAAATTTTTGTTAAATCAGCTCATTTTTTTAAACCAATAGGCCGAAAATCGGCAAAAATCCCTTATAAAT
CAAAAGAATAGACCGAGATAGGGTTGAGTGTGTTCCAGTTTGGAAACAAGAGTCCACTATTAAGAACGTGGACTCCAA
CGTCAAAGGGCGAAAACCGTCTATCAGGGCGATGGCCACTACGTGAACCATCACCTAATCAAGTTTTTTGGGGTCTG
AGGTGCCGTAAAGCACTAAATCGGAACCTAAAGGGAGCCCCGATTTAGAGCTTGACGGGGAAAGCCGGCGAACGTG
GCGAGAAAGGAAGGAAGAAAGCGAAAGGAGCGGGCGTAGGGCGCTGGCAAGTGTAGCGGTACGCTGCGCGTAACC
ACCACACCCGCCGCGCTTAATGCGCCGCTACAGGGCGCGTC

Quality Control Plasmid Transfection:

Lipofectamine LTX Transfection in HT1080 Cells
7-16-14 @ 24 hrs 0.7×10^6 c/well in 6 well plate



G1061 pFBAAVTREmiSafeBgHpA-
CMVrtTAwtIRESeGFP-SV40pA
2ug/well



Quality Control Plasmid ITR Digest:

7-10-14 miniprep 1 @ 432ng/ul

Expected fragments:

1. 1Kb+ ladder
2. uncut
3. BssHII: 4600, 3586, 269, 85, 74
4. MscI: 3503, 2791, 1288, 1032
5. XmaI: 4652, 2256, 1695, 11

Aldevron Maxi Lot 51905 7/17/14



Background:

This vector was designed with two independent cassettes. 1. TRE-miRNA-BgHpA and 2. CMV-rtTA-wtIRESegFP. The dox inducible, tetracycline-regulated TRE promoter drives the miRNA cassette. Constitutive expression from the CMV promoter drives the rtTA transactivator and the expression of eGFP for tracking. In the absence of induction according to the originator Clonetech, the Tet-On Advanced transactivator shows virtually no residual binding to the TRE in P_{Tight} . Thus, basal expression is extremely low.

Please refer to the original Clonetech guides for useful information on working with the Tet-On Advanced System:


Tet-On Advanced Inducible Gene Expression Systems User Manual #PT3898-1_102312-7

http://www.google.com/url?sa=t&rct=j&q=&esrc=s&source=web&cd=7&ved=0CGcQFjAG&url=http%3A%2F%2Fwww.clontech.com%2Fxxclt_ibcGetAttachment.jsp%3FcltemId%3D17561&ei=3LupU6jCJsWfqAbkslLoCA&usq=AFQjCNEQOJMoyr4OPqXaaOPcwNbudctcA&bvm=bv.69620078.d.b2k

Mir-XTM Inducible miRNA Systems User Manual #PT5050-1

http://www.google.com/url?sa=t&rct=j&q=&esrc=s&source=web&cd=1&ved=0CCIQFjAA&url=http%3A%2F%2Fwww.clontech.com%2Fxxclt_ibcGetAttachment.jsp%3FcltemId%3D17566&ei=rMCPu8ftB9KXqAb4mYGAAg&usq=AFQjCNGq-fjnWOPOihXdiOQ0VTfaHbdVw&bvm=bv.69620078.d.b2k

Hairpin Structure:



mU6miSafe cassette EcoRI mU6 promoter miSafe Pol III Terminator EcoRI

GAATTCGAGCGCCGCATCTCTAGGCCCGCGCCGCCCCCTCGCACAGACTTGTGGGAGAAGCTCGGCTACTCCCCTGCCCGGTTAATTTGCATATAATA
 TTCTAGTAAGTATAGAGGCTTAATGTGCGATAAAAGACAGATAATCTGTTCTTTTAAATACTAGCTACATTTACATGATAGGCTTGGATTCTA
 TAAGAGATACAAATACTAAATATTATTTAAAAACAGCACAAAAGAACTCACCTAACTGTAAGTAATTGTGTGTTTGGACTATAAATAT
 CCCTTCGAGAAAAGCCTTGTTCGTTTAGTGAACCGTCAGATGGTACCGTTTAAACTCGAGTGAGCGCAGCGAACGACTTACGCGTTTACTGTAAAG
 CCACAGATGGGTAAACGCGTAAGTCGTTTCGCTACGCCTACTAGA GCGGCCGCCACAGCGGGGAGATCCAGACATGATAAGATACTTTTTGAATTC

miSafe Hairpin

XhoI linker antisense common loop sense linker SpeI/XbaI
 CTCGAGTGAGCGCAGCGAACGACTTACGCGTTTACTGTAAAGCCACAGATGGGTAAACGCGTAAGTCGTTTCGCTACGCCTACTAGA

miSafe designed by Ryan Boudreau and Beverly L. Davidson

aaaactcgagtgagcgcagcgaacgacttacgcggttactgtaaagccacagatgggtaaacgcgtaagtcgctacgccc
 tactagttttt

	A	C		GUAAAG
5'...	AGUG	GCG	AGCGAACGACUUACGCGUUUACU	C
3'...	UCAU	CGC	UCGCUUGCUGAAUGCGCAAUGG	C
	C	A		GUAGACA

Oligo 1: aaaactcgagtgagcgcagcgaacgacttacgcggttactgtaaagccacagatggg

Oligo 2: aaaaactagtaggcgtagcgaacgacttacgcggttaccatctgtggccttacag
 "ctgtaaagccacagatgggtaaacgcgtaagtcgctacgcctactagttttt"

Cumulative Distribution Antisense Seed Complements: TACGCGTT, ACGCGTTT, CGCGTTTA

Information and Cloning Suggestions for Working with pFBAAV Plasmids:

Characteristics of AAV:

Adeno-Associated virus (AAV) is a non-pathogenic helper dependant parvovirus. This is one of the most promising vehicles for gene delivery. Recombinant AAV vectors have predominantly episomal gene expression. AAV has long term expression in terminally differentiated cells. It can infect a wide range of cells, including non-dividing cells.

Insert Size:

The major disadvantage of the AAV vectors is the small genome, which limits the size of the transgene (from ITR to ITR) to about 4.7 Kb. Please consider the size of the promoter, polyA, any reporters, and your gene of interest when planning the cloning into an AAV shuttle vector. When the length of inserted DNA between the 2 ITRs is close the maximal allowed, i.e., 4.7Kb, the packaging efficiency decreases significantly.

AAV Production:

The University of Iowa Viral Vector Core uses two systems for AAV production: A Triple Transfection System and the Baculovirus System.

The Triple Transfection System may be used with any pAAV plasmid including the pFBAAV plasmids. A pHelper plasmid isolated from adenovirus, a Rep/Cap plasmid, and the AAV ITR plasmid package containing the gene of interest are transfected into 40 x 150mm plates of adherent 293FT cells per prep. This system is more labor intensive for production and purification but can be completed within three to four weeks of receiving the plasmid.

The services provided for production are:

- Transfection.
- Large-scale amplification of AAV vectors.
- The AAV vector is purified using either CsCl or an iodixanol gradient followed by ion-exchange or affinity chromatography.
- A physical titer in viral genomes per ml (vg/ml) is assessed by QPCR.
- **Sample required:** 1500ug of high quality AAV plasmid expressing the gene of interest at a concentration greater than 0.25ug/ul per prep. One disadvantage to the investigator in the triple transfection system is having to provide 1500ug of DNA for each subsequent prep and the possibility of recombination during the maxi prep process. Each maxi prep should be digested and evaluated for integrity.

The Baculovirus System can only be used with pFBAAV plasmids with Tn7 transposable elements to incorporate the ITRs and gene of interest into the Bacmid DNA. *Please note: not all AAV plasmids contain the Baculo elements so this system may not be an option for everyone. This system also relies on the gentamicin selection marker within the transposable elements. The system was developed by Dr. Robert Kotin at the National Institute of Health. The system is described in the following publication: "*Insect Cells as a Factory to Produce Adeno-Associated Virus Type 2 Vectors*" [Human Gene Therapy 13:1935-1943](#) (November 1, 2002). The Baculovirus System takes 5-6 weeks for production but has many advantages. The Baculovirus System can tolerate more toxic proteins driven by mammalian promoters as these promoters will not be active in the insect cells. This system is also generally more robust with higher titers, cost-effective, and less labor intensive for production. Once we have a Baculovirus stock of the investigator package, we can make additional preps readily and we have tested the Baculovirus stock to be stable for years.

The services provided for production are:

- Transposition of the gene interest into the DH10Bac bacmid DNA. Everything in the pFBAAV plasmid between the transposable elements (Tn7, ITR, Promoter, gene of interest, pA, ITR, gentamicin resistance, Tn7) is incorporated into the large 135,000kb bacmid DNA. Your plasmid must be gentamicin resistant for the transposition to work.
- Growth and midi prep preparation of Bacmid DNA.
- Transfection of the Bacmid DNA with gene of interest into Sf9 cells. Approximately 12mls of the P1 stage Baculovirus is harvested 5 days later.
- Titer of the Baculovirus P1 stock. We are able to make 1-3 preps of final AAV virus using the P1 stock depending on the titer of the stock and the desired MOI of infection for the final AAV prep.

- For multiple AAV virus preps the P1 stock is amplified to a larger volume P2 Baculovirus stock and titered.
- The final AAV virus is prepared in bioreactor bags on a rocker platform by dual infection with the package Baculovirus and the desired Rep/Cap Baculovirus into Sf9 cells. The Baculovirus provides pHelper function. The culture is harvested 3 days later.
- The AAV vector is purified using either CsCl or an iodixanol gradient followed by ion-exchange or affinity chromatography.
- A physical titer in viral genomes per ml (vg/ml) is assessed by QPCR.
- **Sample required:** 10ug of pFBAAV plasmid expressing the gene of interest at a concentration greater than 0.25ug/ul.

Bacterial Backbone:

All of our internal AAV plasmids contain Baculo compatible elements from Invitrogen pFastBac™. For information on the pFastBac backbone and the Baculovirus System, please review the Invitrogen manual: http://tools.invitrogen.com/content/sfs/manuals/bactobac_topo_exp_system_man.pdf

**Please note: not all AAV plasmids are compatible with the Baculo system.

ITRs (Inverted Terminal Repeats):

The ITRs for AAV were engineered into the Invitrogen pFastBac backbone by Dr. Robert Kotin.

Orientation and Updates:

The MCS and BgHpA were inserted in the antisense direction relative to the Invitrogen backbone and ITRs. Please note that the backbone and ITRs were re-sequenced and updated in 2012 and there are several changes. Please discard any older sequences. The nomenclature of 5'ITR (130bp) and 3'ITR (141bp) have been retained on some maps to identify the sequences. The orientation of the ITRs makes no difference to the formation of the viral vector.

Recombination:

Recombination is a possibility at both the transposable elements (Tn7L and Tn7R) and the Inverted Terminal Repeats (ITRs) in the pFBAAV plasmid. In order to check for and reduce the chance of recombination, we have the following recommendations.

Antibiotic Resistance:

The gentamicin resistant gene lies within the transposable elements and may be lost if the transposable elements recombine. We suggest using **both ampicillin and gentamicin** selection on plates and in your liquid media throughout the cloning process and during amplification to midi or maxi preps.

Ampicillin concentration: 100ug/ml final

Gentamicin Sulfate concentration: 7ug/ml final

E. Coli Competent Cell Recommendations:

We recommend using a stable E. coli strain such as **SURE2, Stable2, or Stable3** for transformation of your final plasmid product. These competent cell strains have been engineered to stop unwanted rearrangement events and lack the components of the pathways that catalyze the rearrangement and deletion of nonstandard secondary and tertiary structures, including cruciforms (caused by inverted repeats) and Z-DNA, that occur frequently in eukaryotic cells. Cloning, however, can be difficult in these strains. We recommend DH5a competent cells for sub-cloning.

Quality Control ITR Digest:

We require a quality control digest of the plasmid upon submission to ensure that there is no recombination before vector production is started. The ITRs in the pFBAAV plasmids have several convenient restriction sites to determine whether the ITRs are intact without sequencing. Sequencing the ITRs can be difficult due to their hairpin secondary structure. We suggest a single digest each of **BssHII, MscI, and XmaI**. It is a good idea to also run un-cut and linearized plasmid next to these digests. Determine your expected fragment sizes for each digest and check carefully to see that you get what you expect and have no extraneous bands. Be sure to check your final midi or maxi product (not just your miniprep) as recombination is possible during any amplification process.

Sequencing pFBAAV plasmids and ITRs

The majority of the plasmid and insert can be sequenced normally. Sequencing of the ITRs is difficult but possible. Sequencing of the ITRs or other difficult secondary structures may be aided by requesting "Extended Denaturation" or by requesting "Betaine Chemistry for AAV ITRs" in the comments section at the University of Iowa Sequencing Core. Non-UI investigators should consult with their own sequencing facility for help if sequencing of the ITRs is desired. We routinely request "Extended Denaturation" for all AAV plasmid sequencing.

There is no one set of primers that will always sequence through the ITRs. The ITR primers listed below have been used with success in many cases. Try sequencing in both directions: from the plasmid in and from your insert out.

Suggested Sequencing Primers for a variety of pFBAAV plasmids:

Below is a general list of commonly used primers at the University of Iowa Viral Vector Core. Please check to make sure the primers work for the shuttle plasmid you have received. There are many free primer design programs available to design primers that will sequence your insert.

BgHpA For	TGAGGAAATTGCATCGCATTGTCT
BgHpA Rev	AGGAAAGGACAGTGGGAGTG
pFBAAV (130bp ITR)	GCCTTGCTGTTCTTCTACGG
pFB587F (141bp ITR)	CTCTACAAATGTGGTATGGCTG
pFB214F (141bp ITR)	GGGGTGGAAATGGAGTTT
pFB158F (141bp ITR)	CATAACAGGAAGAAAAATGCCCCG
EGFP-141Rev	GAACTTCAGGGTCAGCTTGCCGTA
EGFP-651For	TCACATGGTCCTGCTGGAGTT
IRES Rev	CTCACATTGCCAAAAGACG
mU6 For	ACAGACTTGTGGGAGAAGC
CMV Rev	GTAGGAAAGTCCCATAAGGTCA
CMV For	GTGGGAGGTCTATATAAGCAGAGCTCG
IRESfor	GTTGTGAGTTGGATAGTTGTGG

Please contact us with any questions:

Viral Vector Core

vectors@uiowa.edu

University of Iowa

500 Newton Road

221 Eckstein Medical Research Building

Iowa City, IA 52242

Tel: (319) 335-6726

<http://www.medicine.uiowa.edu/vectorcore/>

Updated 2-22-17 AN