

Eyelid

Special arrangements

No special arrangements are necessary.

Fixation

10% neutral buffered formalin.

Special instructions

- ❖ If the resection margins are important, please make a sketch on the Ocular Pathology Consultation Request form as to the location of the tumor or area of interest in the eyelid so that the resection margins can be distinguished in the specimen.
- ❖ Follow the instructions listed under General Information for Routine Specimen Collection for proper packaging and labeling of routine specimens and how to fill out the Ocular Pathology request form.
- ❖ Follow the procedure under General Information on Collection and Handling of Specimens for proper specimen sending information.