

Iris

Special arrangements

If the iris is removed for a suspected tumor, please see instructions listed above under **ciliary body (iridocyclectomy)** specimens.

Fixation

10% neutral buffered formalin.

Special instructions

- ❖ There are no special instructions if the iris is removed at the time of glaucoma filtering procedures.
- ❖ Follow the instructions listed under General Information for Routine Specimen Collection for proper packaging and labeling of routine specimens and how to fill out the Ocular Pathology request form.
- ❖ Follow the procedure under General Information on Collection and Handling of Specimens for proper specimen sending information.